

Anti-D-Lactic acid-D-Lactic acid-Histone H3-K9 antibody (STJ11107135)

STJ11107135

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	DB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	DB:1:500-1:1000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	8350/8351/8352/8353/8354/8355/8356/8357/8358/8968 8290
Gene Symbol	H3C1.H3C2.H3C3.H3C4.H3C6.H3C7.H3C8.H3C10.H3C11.H3C12 H31_HUMAN H31T_HUMAN
Immunogen	ARKST
Immunogen Region	
Specificity	A synthetic D-lactylated peptide around K9 of human Histone H3 (NP_003520.1).
Immunogen Sequence	ARKST