

Anti-NSUN7 antibody (152-280 aa) (STJ11107137)

STJ11107137

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000
Formulation	PBS with 0.09% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	70918
Gene Symbol	Nsun7
Uniprot ID	NSUN7_MOUSE
Immunogen	
Immunogen Region	152-280 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 152-280 of mouse NSUN7 (NP_001387843.1).
Immunogen Sequence	NEETIPEVQEVENLLNGFKT KLAALARCRIKHDALSIYH ILPETVRKQEQRASLTPLYA WINTSKISLEEYNNLRKKG YSKVKSITSVNEKVYAVDQH CFNVLIFPAHLKTDLLNIDL IKDYKLIFQ