

Anti-SMARCB1/SNF5 antibody (286-385 aa) (STJ11107138)

STJ11107138

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	6598
Gene Symbol	SMARCB1
Uniprot ID	SNF5_HUMAN
Immunogen	
Immunogen	286-385 aa
Region	
Specificity	A synthetic peptide corresponding to a sequence within amino acids 286-385 of human SMARCB1/SNF5 (NP_003064.2).
Immunogen	KENSPEKFALKLCSELGLGG EFVTTIAYSIRGQLSWHQKT YAFSENPLPTVEIAIRNTGD ADQWCPLLETLTDAEMEKKI
Sequence	RDQDRNTRRMRLANTAPAW