

Anti-MHC Class II-I-A/I-E antibody [SRM] (STJ11107143)
STJ11107143

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	FC/ELISA
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	FC:1:500-1:1000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	NA NA NA NA NA NA NA
Gene Symbol	H2-Aa H2-Eb1 HA2F_MOUSE HB2I_MOUSE HB2Q_MOUSE
Immunogen	
Immunogen Region	
Specificity	Recombinant protein of mouse MHC Class II (I-A/I-E).
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081