

Anti-SMARCB1/SNF5 antibody (50-150 aa) [SRM] (STJ11107145)

STJ11107145

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:2000-1:10000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	6598
Gene Symbol	SMARCB1
Uniprot ID	SNF5_HUMAN
Immunogen	
Immunogen	50-150 aa
Region	
Specificity	A synthetic peptide corresponding to a sequence within amino acids 50-150 of human SMARCB1/SNF5 (NP_003064.2).
Immunogen	LWRRLATVEERKKIVASSHG KKT KPNTKDHGYTTLATSVT LLKASEVEEILDGNDEKYKA VSISTEPPTYLREQKAKRNS
Sequence	QWVPTLPNSSHHLDVPCST T