

Anti-KLF3 antibody (186-256 aa) (STJ11107146)

STJ11107146

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	16599
Gene Symbol	KLF3
Uniprot ID	KLF3_MOUSE
Immunogen	
Immunogen Region	186-256 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 186-256 of mouse KLF3 (NP_032479.1).
Immunogen Sequence	SYEKPLLQKKIKIEPGIEPQ RTDYYPEEMSPPLMNPVSPP QALLQENHPSVIVQPGKRPL PVESPDTRKR