

Anti-ZDHHC9 antibody (254-364 aa) (STJ11107147)

STJ11107147

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	208884
Gene Symbol	Zdhhc9
Uniprot ID	ZDHC9_MOUSE
Immunogen	
Immunogen Region	254-364 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 254-364 of mouse ZDHHC9 (NP_766053.1).
Immunogen Sequence	QTTNEDIKGSWTGKNRVQNP YSHGNIVKNCCEVLCGPLPP SVLDRRGILPLEESGSRPPS TQETSSLLPQSPASTEHMN SNEMAEDTSIPEEMPPPEPP EPPQEASEAEK