

Anti-NF Kappa B P105 antibody (741-968 aa) [SRM] (STJ11107156)

STJ11107156

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/IP/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	IP:0.5 Mu g-4 Mu g antibody for 200 Mu g-400 Mu g extracts of whole cells ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	4790
Gene Symbol	NFKB1
Uniprot ID	NFKB1_HUMAN
Immunogen	
Immunogen Region	741-968 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 741-968 of human NF Kappa B P105 (NP_001158884.1).
Immunogen Sequence	AGADPLVENFEPLYDLDDSW ENAGEDEGVVPGTTPLDMAT SWQVFDILNGKPYEPEFTSD DLLAQGDMKQLAEDVKLQLY KLLLEIPDPDKNWATLAQKLG LGILNNAFRLSPAPSKTLM D NYEVSGGTVRELVEALRQMG YTEAIEVQAASSPVKTT SQ AHSLLPSASTRQQIDELRD SDSVCDSGVETSFRKLSFTE SLTSGASLLTLNKMPhDYGQ EGPLEGKI