

## Anti-ACTN1 antibody (281-380 aa) [SRM] (STJ11107157)

STJ11107157

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	
<b>Applications</b>	WB/ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human/Mouse

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	SRM
<b>Concentration</b>	Lot specific
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity purification
<b>Dilution Range</b>	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
<b>Formulation</b>	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
<b>Isotype</b>	IgG
<b>Storage</b>	
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">87</a>
<b>Gene Symbol</b>	<a href="#">ACTN1</a>
<b>Uniprot ID</b>	<a href="#">ACTN1_HUMAN</a>
<b>Immunogen</b>	
<b>Immunogen</b>	281-380 aa
<b>Region</b>	
<b>Specificity</b>	A synthetic peptide corresponding to a sequence within amino acids 281-380 of human Alpha-Actinin (NP_001093.1 ).
<b>Immunogen</b>	KLASDLLEWIRRTIPWLENR VPENTMHAMQQKLEDFRDYR RLHKPPKVQEKCQLEINFNT LQTKLRLSNRPAFMPSEGRM
<b>Sequence</b>	VSDINNAWGCLEQVEKGYEE