

Anti-Aurora A antibody (1-100 aa) [SRM] (STJ11107162)

STJ11107162

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	6790
Gene Symbol	AURKA
Uniprot ID	AURKA_HUMAN
Immunogen	
Immunogen Region	1-100 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Aurora A (NP_003591.2).
Immunogen Sequence	MDRSKENCISGPVKATAPVG GPKRVLVTQQFPCQNPLPVN SGQAQRVLCPSNSSQRIPLQ AQKLVSSHKPVQNQKQKQLQ ATSVPHPVSRPLNNTQKSKQ