

Anti-CDC20 antibody (1-100 aa) [SRM] (STJ11107164)
STJ11107164**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	
Applications	WB/IP/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 IP:0.5 Mu g-4 Mu g antibody for 200 Mu g-400 Mu g extracts of whole cells ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	991
Gene Symbol	CDC20
Uniprot ID	CDC20_HUMAN
Immunogen	
Immunogen Region	1-100 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CDC20 (NP_001246.2).
Immunogen Sequence	MAQFAFESDLHSLQLDAPI PNAPPARWQRKAKEAAGPAP SPMRAANRSHSAGRTPGRTP GKSSSKVQTTPSKPGGDRYI PHRSAAQMEVASFLLSKENQ