

Anti-COP1/RFWD2 antibody (582-731 aa) [SRM] (STJ11107166)

STJ11107166

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	IHC-P:1:50-1:200
	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	64326
Gene Symbol	COP1
Uniprot ID	COP1_HUMAN
Immunogen	
Immunogen Region	582-731 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 582-731 of human COP1/RFWD2 (NP_071902.2).
Immunogen Sequence	YYDLRNTKQIMVFKGHRKA VSYAKFVSGEEIVSASTDSQ LKLWNVGKPYCLRSFKGHIN EKNFVGLASNGDYIACGSEN NSLYLYYKGLSKTLLTFKFD TVKSVLDKDRKEDDTNEFVS AVCWRALPDGESNVLIAANS QGTIKVLELV