

## Anti-GABPB1 antibody (290-395 aa) [SRM] (STJ11107175)

STJ11107175

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	
<b>Applications</b>	WB/ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	SRM
<b>Concentration</b>	Lot specific
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity purification
<b>Dilution Range</b>	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
<b>Formulation</b>	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
<b>Isotype</b>	IgG
<b>Storage</b>	
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	2553
<b>Gene Symbol</b>	GABPB1
<b>Uniprot ID</b>	GABP1_HUMAN
<b>Immunogen</b>	
<b>Immunogen Region</b>	290-395 aa
<b>Specificity</b>	Recombinant fusion protein containing a sequence corresponding to amino acids 290-395 of human GABPB1 (NP_005245.2).
<b>Immunogen Sequence</b>	TSGIGQPIIVTMPDGQQVLT VPATDIAEETVISEEPPAKR QCIEIENRVESAEIEEREA LQKQLDEANREAQKYRQQLL KKEQEAAYRQKLEAMTRLQ TNKEAV