

Anti-Grp94 antibody (616-715 aa) [SRM] (STJ11107176)
STJ11107176**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	7184
Gene Symbol	HSP90B1
Uniprot ID	ENPL_HUMAN
Immunogen	
Immunogen Region	616-715 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 616-715 of human Grp94 (NP_003290.1).
Immunogen Sequence	EPLLNWMKDKALKDKIEKAV VSQRLTESPCALVASQYGWS GNMERIMKAQAYQTGKDIST NYYASQKKTFEINPRHPLIR DMLRRIKEDEDDKTVLDLAV