

Anti-EPOR antibody (25-249 aa) [SRM] (STJ11107179)

STJ11107179

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	2057
Gene Symbol	EPOR
Uniprot ID	EPOR_HUMAN
Immunogen	
Immunogen Region	25-249 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 25-249 of human EPOR (NP_000112.1).
Immunogen Sequence	APPPNLPDPKFESKAALLAA RGPEELLCFTERLEDLVCFW EEASAGVGPGNYSFSYQLE DEPWKLCRLHQAPTARGAVR FWCSLPTADTSSFVPLELRV TAASGAPRYHRVIHINEVVL LDAPVGLVARLADESGHVVLRWLPPPETPMTSHIRYEVDV SAGNGAGSVQRVEILEGRTE CVLSNLRGRTRYTFVRARM AEPSFGGFWSAWSEPVSLLT PSDLD