

Anti-TNFRSF11A antibody (30-212 aa) [SRM] (STJ11107180)

STJ11107180

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	IF/ICC:1:50-1:200
	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	8792
Gene Symbol	TNFRSF11A
Uniprot ID	TNR11_HUMAN
Immunogen	
Immunogen Region	30-212 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 30-212 of human TNFRSF11A (NP_003830.1).
Immunogen Sequence	IAPPCTSEKHYEHLGRCCNK CEPGKYMSSKCTTSDSVCL PCGPDEYLDWNEEDKCLLH KVCDTGKALVAVVAGNSTTP RRCACTAGYHWSQDCECCRR NTECAPGLGAQHPLQLNKDT VCKPCLAGYFSDAFSSTDKC RPWTNCTFLGKRVEHHGTEK SDAVCSSSLPARKPPNEPHV YLP