

Anti-INTS12 antibody (1-230 aa) [SRM] (STJ11107182)

STJ11107182

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	57117
Gene Symbol	INTS12
Uniprot ID	INT12_HUMAN
Immunogen	
Immunogen Region	1-230 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-230 of human INTS12 (NP_065128.2).
Immunogen Sequence	MAATVNLELDPIFLKALGFL HSKSKDSAEKLKALLDESLA RGIDSSYRPSQKDVEPPKIS STKNISIKQEPKISSSLPSG NNNGKVLTTKVKKEAEKRP ADKMKSDITEGVDPKPKPRL EKPETQSSPITVQSSKDLPM ADLSSFEETSADDFAMEMGL ACVVCQRQMMVASGNQLVECQ ECHNLYHRDCHKPQVTDKEA NDPRLVWYCARCTRQMKRMA QKTQKPPQKP