

Anti-LRPPRC antibody (1014-1394 aa) [SRM] (STJ11107183)

STJ11107183

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	IHC-P:1:50-1:200 IF/ICC:1:50-1:200 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	10128
Gene Symbol	LRPPRC
Uniprot ID	LRPPRC_HUMAN
Immunogen	
Immunogen Region	1014-1394 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1014-1394 of human LRPPRC (NP_573566.2).
Immunogen Sequence	LWYEDEKHSLSNSTTTP DFQKDILJACRLNQKKGAYD IFLNAKEQNIVFNAETYSNL IKLLMSSEDYFTQAMEVKAFAL ETHIKGFTLNDAANSRLIIT QVRRDYLKEAVTTLKTVLDQ QQTPSRLAVTRVIQALAMKG DVENIEVVQKMLNGLEDSIG LSKMVFINNIALAQIKNNNI DAAIENIENMLTSENKVIEP QYFGLAYLFRKVIEEQLEPA VEKISIMAERLANQFAIYK