

Anti-LTBR antibody (101-200 aa) [SRM] (STJ11107184)

STJ11107184

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	4055
Gene Symbol	LTBR
Uniprot ID	TNR3_HUMAN
Immunogen	
Immunogen	101-200 aa
Region	
Specificity	A synthetic peptide corresponding to a sequence within amino acids 101-200 of human LTBR (NP_002333.1).
Immunogen	CRPCDPVMGLEEIAPCTSKR KTQCRCQPGMFCAAWALECT HCELLSDCPPGTEAEKDEV GKGNNHCVPCKAGHFQNTSS
Sequence	PSARCQPHTRCENQGLVEAA