

Anti-CCL11/Eotaxin antibody (26-97 aa) [SRM] (STJ11107189)

STJ11107189

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	20292
Gene Symbol	Ccl11
Uniprot ID	CCL11_MOUSE
Immunogen	
Immunogen Region	26-97 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 26-97 of mouse CCL11/Eotaxin (NP_035460.1).
Immunogen Sequence	GSIPTSCCFIMTSKKIPNTL LKSYKRITNNRCTLKAIVFK TRLGKEICADPKKKVWQDAT KHLQKLQTPKP