

Anti-MRC2 antibody (35-513 aa) [SRM] (STJ11107192)

STJ11107192

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:1000-1:2000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	9902
Gene Symbol	MRC2
Uniprot ID	MRC2_HUMAN
Immunogen	
Immunogen Region	35-513 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 35-513 of human MRC2 (NP_006030.2). DAALPEPNVFLIFSHGLQGC LEAQGGQVRVTPACNTSLPA QRWKWVSRNRLFNLGTMQCL GTGWPGTNTTASLGMYESCDR EALNLRWHCRTLGDQLSLLL GARTSNISKPGTLERGDQTR SGQWRIYGSEEDLCALPYHE VYTIQGNSHGKPCPTIPFKYD NQWFHGCSTSTGREDDHLWCA TTQDYGKDERWGFCPIKSND CETFWDKDLTDSCYQFNQ STLSWREAWASCEQQGADL



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081