

Anti-SCN3B antibody (25-145 aa) [SRM] (STJ11107201)

STJ11107201

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 IHC-P:1:50-1:200 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	55800
Gene Symbol	SCN3B
Uniprot ID	SCN3B_HUMAN
Immunogen	
Immunogen Region	25-145 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 25-145 of human SCN3B (NP_001035241.1).
Immunogen Sequence	VCVEVPSETEAVQGNPMKLR CISCMKREEVEATTVVEWYF RPEGGKDFLIYEYRNGHQEV ESPFQGRQLQWNGSKDLQDVS ITVLNVTLNDSGLYTCNVSREFEFEAHRPFVKTTRLIPLR V