

Anti-Syntaxin 6 antibody (1-234 aa) [SRM] (STJ11107202)

STJ11107202

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	10228
Gene Symbol	STX6
Uniprot ID	STX6_HUMAN
Immunogen	
Immunogen Region	1-234 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-234 of human Syntaxin 6 (NP_001273139.1).
Immunogen Sequence	MSMEDPFFVVKGEVQKAVNT AQGLFORWTELLQDPSTATR EEIDWTTNELRNLRNLSIEWD LEDLDETISIVEANPRKFNLDATLSIRKAFITSTRQVVR DMKDDQMSSTSSVQALAEKRNQ QALLGDSGSQNWSTGTTDKY GRLDRELQRANSHFIEEQQAQQQLIVEQQDEQLELVSGSI GVLKNMSQRIGGELEEQAVM LEDFSHELESTQSRLDNVMK KLAKVSHMTSDRRQ