

Anti-TAF1C antibody (75-221 aa) [SRM] (STJ11107203)

STJ11107203

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	9013
Gene Symbol	TAF1C
Uniprot ID	TAF1C_HUMAN
Immunogen	
Immunogen Region	75-221 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 75-221 of human TAF1C (NP_005670.4).
Immunogen Sequence	PWDPGLTARDLLFRGGCRYR KRPRVLDVTEQISRFLDH GDVAFAPLGKLMLENFKLEG AGSRTKKKT VSVKKLLQDL
Sequence	GGHQPWGCPWAYLSNRQRRF SILGGPILGTSVASHLAELL HEELVLRWEQLLLDEACTGG ALAWVPG