

Anti-ZIP14 antibody (31-157 aa) [SRM] (STJ11107205)
STJ11107205**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 IHC-P:1:50-1:200 IF/ICC:1:50-1:200 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	23516
Gene Symbol	SLC39A14
Uniprot ID	S39AE_HUMAN
Immunogen	
Immunogen Region	31-157 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 31-157 of human ZIP14 (NP_001121903.1).
Immunogen Sequence	SSLGAPAI SAASFLQDLIHR YGEGDSLTLQQLKALLNHL D VGVGRGNVTQH VQGHRLST CFSSGDLFTA HNFSEQSRIG SSELQEF CPTILQQ LDSRAC TSENQENEENEQTEEGRPSA VEWGYG