

Anti-CD163 antibody [SRM] (STJ11107207)

STJ11107207

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:1000-1:9000 IHC-P:1:10000-1:40000 IF/ICC:1:2000-1:4000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	93671 9332
Gene Symbol	Cd163 CD163
Uniprot ID	C163A_MOUSE C163A_HUMAN
Immunogen	
Immunogen Region	
Specificity	Recombinant protein of human/mouse CD163.
Immunogen Sequence	