

Anti-CD66a/CEACAM1 antibody [SRM] (STJ11107209)
STJ11107209**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:3000-1:12000 IHC-P:1:4000-1:10000 IF/ICC:1:2000-1:8000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	26365 634
Gene Symbol	Ceacam1 CEACAM1
Uniprot ID	CEAM1_MOUSE CEAM1_HUMAN
Immunogen	
Immunogen Region	
Specificity	Recombinant protein of human/mouse CD66a/CEACAM1.
Immunogen Sequence	