

Anti-Ki67 antibody (1816-2163 aa) [SRM] (STJ11107215)

STJ11107215

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:1000-1:5000
Range	IHC-P:1:2000-1:12000 IF/ICC:1:2000-1:5000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	17345
Gene Symbol	Mki67
Uniprot ID	Ki67_MOUSE
Immunogen	
Immunogen Region	1816-2163 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1816-2163 of mouse Ki67 (NP_001074586.2).
Immunogen Sequence	DETTKIALQSPQPGHIINPA SMKRQSNMSLRKDMREFSIL EKQTQSRGRDAGTPAPMQEE NGTTAIMETPKQKLDFIGNS TGHKRRPRTPKNRAQPLEDL DGFQELFQTPAGASDPVSVE ESAKISLASSQAEPV RTPAS TKRRSKTGLSKVDVRQEPST LGKRMKSLGRAPGTPAPVQE ENDSTAFMETPKQKLDFTGN SSGHKRRPQTPKIRAQPLED LDGFQELFQTPAGANDSVT