

Anti-Granulin antibody (494-593 aa) [SRM] (STJ11107218)

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/IF/CC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:1000-1:3000 IF/CC:1:200-1:800 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	2896
Gene Symbol	GRN
Uniprot ID	GRN_HUMAN
Immunogen	
Immunogen Region	494-593 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 494-593 of human Granulin (NP_002078.1).
Immunogen Sequence	SCEKEVWSAQPATFLARSPH VGVKDVECGEGHFCHDNQTC CRDNRQGWACCPYRQGVCCA DRRHCCPAGFRCAARGTKCL RREAPRWDAPLRDPALRQLL