

Anti-TAZ antibody (1-270 aa) [SRM] (STJ11107221)

STJ11107221

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:1500-1:9000
Range	IHC-P:1:50-1:200
Formulation	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Isotype	PBS with Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Storage	IgG
Instruction	

TARGET INFORMATION

Gene ID	25937
Gene Symbol	WWTR1
Uniprot ID	WWTR1_HUMAN
Immunogen	
Immunogen Region	1-270 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human TAZ (NP_056287.1)
Immunogen Sequence	MNPASAPPLPPPGQVIHV TQDLDTLEALFNSVMNPKP SSWRKKILPESFFKEPDGSGS HSRQSSTDSSGGHPGPRLAG GAQHVRSHSSPASLQLGTGA GAAGSPAQQHAHLRQQSYDV TDELPLPPGWEMTFTATGQR YFLNHIEKITTWQDPRKAMN QPLNHMNLHPAVSSTPVPQR SMAVSQPNLVMNHQHQQQMA PSTLSQQNHPTQNPAGLMS MPNALTTQQQQQKLRLQR