

Anti-GLDC antibody (580-749 aa) [SRM] (STJ11107222)

STJ11107222

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:25000-1:150000
Range	IHC-P:1:4000-1:16000 IF/ICC:1:2000-1:8000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	2731
Gene Symbol	GLDC
Uniprot ID	GCSF_HUMAN
Immunogen	
Immunogen Region	580-749 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 580-749 of human GLDC (NP_000161.2).
Immunogen Sequence	HPFVPLDQAQGYQQLFRELE KDLCELTGYDQVCFQPNSSA QGEYAGLATIRAYLNQKGEG HRTVCLIPKSAHGTNPASAH MAGMKIQPVEVDKYGNIDAV HLKAMVDKHKENLAIMITY PSTNGVFEENISDVCDLIHQ HGGQVYLDGANMNAQVGICR PGDFGSDVSH