

Anti-CD69 antibody (62-199 aa) [SRM] (STJ11107231)

STJ11107231

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:2500-1:10000
Range	IHC-P:1:200-1:800 IF/ICC:1:200-1:800 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	969
Gene Symbol	CD69
Uniprot ID	CD69_HUMAN
Immunogen	
Immunogen Region	62-199 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 62-199 of human CD69 (NP_001772.1).
Immunogen Sequence	SVGQYNCPGQYTFSPSDSH VSSCEDWVG YQRKCYFIST VKRSWTS AQNACSEHGATLA VIDSEKDMNFLKRYAGREEH WVGLKKEPGHPWKWSNGKEF NNWFNVTGSDKCVFLKNTEV SSMECEKNLYWICNKPYPK