

Anti-UBA1/7 antibody (959-1058 aa) [SRM] (STJ11107241)

STJ11107241

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:1000-1:10000 IHC-P:1:500-1:2000 IF/ICC:1:200-1:800 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	7318 7317
Gene Symbol	UBA7 UBA1
Uniprot ID	UBA7_HUMAN UBA1_HUMAN
Immunogen	
Immunogen Region	959-1058 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 959-1058 of human UBA1/7 (P22314).
Immunogen Sequence	EVQGLQPNGEEMTLKQFLDY FKTEHKLEITMLSQGVSMY SFFMPAAKLKERLDQPMTEI VSRVSKRKLGRHVRALVLEL CCNDESGEDVEVPYVRYTIR