

Anti-CDK1/2 antibody (199-298 aa) [SRM] (STJ11107242)

STJ11107242

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:3000-1:12000 IHC-P:1:200-1:2000 IF/ICC:1:200-1:800 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	983 1017
Gene Symbol	CDK1 CDK2
Uniprot ID	CDK1_HUMAN CDK2_HUMAN
Immunogen	
Immunogen Region	199-298 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 199-298 of human CDK1/2 (P24941).
Immunogen Sequence	RRALFPGDSEIDQLFRIFRT LGTPDEVVWPGVTSMPTYKP SFPKWARQDFSKVVPPLDED GRSLLSQMLHYDPNKRISAK AALAHPPFFQDVTKPVPHLRL