

Anti-CD11b antibody (1-100 aa) [SRM] (STJ11107243)
STJ11107243**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:1000-1:6000 IHC-P:1:200-1:800 IF/ICC:1:200-1:800 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	3684
	NA
Gene Symbol	ITGAM
	Itgam
Uniprot ID	ITAM_HUMAN
	ITAM_MOUSE
Immunogen	
Immunogen	1-100 aa
Region	
Specificity	Recombinantfusion protein containing a sequence corresponding to amino acids 1-100 of human CD11b (NP_000623.2).
Immunogen	MALRVLLLTALTLCGHNLD TENAMTFQENARGFGQSVVQ LQGSRRVVVGAPQEIVAANQR GSLYQCDYSTGSCEPIRLQV
Sequence	PVEAVNMSLGLSLAATTSP