

Anti-Ki67 antibody (2271-2370 aa) [SRM] (STJ11107244)

STJ11107244

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:1000-1:6000 IHC-P:1:1000-1:4000 IF/ICC:1:200-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	4288 17345
Gene Symbol	MKI67 Mki67
Uniprot ID	KI67_HUMAN KI67_MOUSE
Immunogen	
Immunogen Region	2271-2370 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 2271-2370 of human Ki67 (NP_002408.3).
Immunogen Sequence	PAGGDEKDMKAFMGTPVQKL DLPGNLPGSKRWPQTPKEKA QALDLAGFKELFQTPGTDK PTTDEKTTKIACKSPQDPDV DTPASTKQRPKRNLRKADVE