

Anti-SFTPC antibody (91-190 aa) [SRM] (STJ11107253)

STJ11107253

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:10000-1:100000 IHC-P:1:4000-1:16000 IF/ICC:1:200-1:800 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	6440
Gene Symbol	SFTPC
Uniprot ID	PSPC_HUMAN
Immunogen	
Immunogen Region	91-190 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 91-190 of human SFTPC (NP_003009.2).
Immunogen Sequence	TATFSIGSTGLVVYDYQQL IAYKPAPGTCCYIMKIAPES IPSLEALTRKVHNFQMECSL QAKPAVPTSKLGQAEGRDAG SAPSGGDP AFLGMAVSTLCG