

Anti-IRP2 antibody (700-963 aa) [SRM] (STJ11107276)

STJ11107276

GENERAL INFORMATION



Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:2500-1:10000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	3658
Gene Symbol	IREB2
Uniprot ID	IREB2_HUMAN
Immunogen	
Immunogen Region	700-963 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 700-963 of human IRP2 (NP_004127.2).
Immunogen Sequence	SLEAPDSVLFPWDLKSTYIR CPSFFDKLTKEPIALQAIEN AHVLLYLGDSVTTDHISPAQ SIARNSAAAKYLTNRGLTPR EFNSYGARRGNDVMTRGTF ANIKLFNFKFIGKPAPKTIHF PSGQTLDDVFEEALYQKEGI PLIILAGKKYSGNSRDWAA KGPYLLGVKAVLAESYEKIH KDHLIGIGIAPLQFLPGENA DSLGLSGRETFSLTPEELS PGITLNIQTSTGKVFVIA



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081