

Anti-IL33 antibody (109-270 aa) [SRM] (STJ11107280)

STJ11107280

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	IF/ICC:1:200-1:800
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	90865
Gene Symbol	IL33
Uniprot ID	IL33_HUMAN
Immunogen	
Immunogen Region	109-270 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 109-270 of human IL33 (NP_001300974.1).
Immunogen Sequence	HDSSITGISPITEYLASLST YNDQSITFALEDSESYEIVVE DLKKDEKKDKVLLSYYESQH PSNESGDGVDGKMLMVTLSPTKDFWLHANNKEHSVELHKC EKPLPDQAFFVLHNMHSNCV SFECKTDPGVFIGVDNHLA LIKVDSSSENLC TENILFKLS ET