

Anti-PSMC4 antibody (1-165 aa) [SRM] (STJ11107292)

STJ11107292

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	IHC-P:1:200-1:800
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	5704
Gene Symbol	PSMC4
Uniprot ID	PRS6B_HUMAN
Immunogen	
Immunogen Region	1-165 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-165 of human PSMC4 (NP_006494.1).
Immunogen Sequence	MEEIGILVEKAQDEIPALSV SRPQTGLSFLGPEPEDLEDL YSRYKKLQQELEFLEVQEEY IKDEQKNLKKKEFLHAQEEVK
Sequence	RIQSIPLVIGQFLEAVDQNT AIVGSTTGSNYVVRILSTID RELKPNASVALHKHSNALV DVLPEADSSIMMLTSDQKP DVMYA