

Anti-NGLY1 antibody (1-300 aa) [SRM] (STJ11107293)

STJ11107293

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	IHC-P:1:100-1:500
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	55768
Gene Symbol	NGLY1
Uniprot ID	NGLY1_HUMAN
Immunogen	
Immunogen Region	1-300 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human NGLY1 (NP_060767.2).
Immunogen Sequence	MAAAALGSSSGSASPAVAEL CQNTPETFLEASKLLTYAD NILRNPNDEKYRSIRIGNTA FSTRLLPVRGAVECLFEMGF EEGETHLIFPKKASVEQLQK IRDLIAIERSRLDGSNKSH KVKSSQQAASTQLPTTPSS NPSGLNQHTRNRQQGSSDPP SASTVAADSAILVLSQNIQ HVLVYENPALQEALACIPV QELKRKSQEKLRSARKLDKG INISDEDFLLLELLHWFKE



Immunohistochemistry analysis of paraffin-embedded Human kidney tissue using NGLY1 Rabbit mAb (STJ11107293) at a dilution of 1:500 (40x lens). High pressure antigen retrieval performed with 0.01M Citrate Buffer (pH 6.0) prior to immunohistochemistry staining.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081