

Anti-NF2/Merlin antibody (347-595 aa) [SRM] (STJ11107299)

STJ11107299

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:3000-1:12000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	4771
Gene Symbol	NF2
Uniprot ID	MERL_HUMAN
Immunogen	
Immunogen Region	347-595 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 347-595 of human NF2/Merlin (NP_000259.1).
Immunogen Sequence	EEAERTRDELERRLLQMKEE ATMANEALMRSEETADLLAE KAQITEEEAKLLAQKAAEAE QEMQRIKATAIRTEEEKRLM EQKVLEAEVLALKMAEESER RAKEADQLKQDLQEAREAE RAKQKLEIATKPTYPPMNP IPAPLPDPSPFNIGDLSL DFKDTDMKRLSMEIEKEV EYMEKSKHLQEQLNELKTEI EALKLKERETALDILHNENS DRGGSSSKHNTIKKLTLSA