

Anti-CHAT antibody (461-560 aa) [SRM] (STJ11107317)

STJ11107317

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IF/ICC:1:200-1:800 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	1103
Gene Symbol	CHAT
Uniprot ID	CLAT_HUMAN
Immunogen	
Immunogen Region	461-560 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 461-560 of human CHAT (NP_065574.4).
Immunogen Sequence	VTQSSRKLIRADSVSELPAP RRLRWKCSPEIQGHLASSAE KLQRIVKNLDFIVYKFDNYG KTFIKKQKCSFDAFIQVALQ LAFYRLHRRLVPTYESASIR