

Anti-CD62E/E-Selectin antibody (22-230 aa) [SRM] (STJ11107318)

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:1000-1:4000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	6401
Gene Symbol	SELE
Uniprot ID	LYAM2_HUMAN
Immunogen	
Immunogen Region	22-230 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 22-230 of human CD62E/E-Selectin (NP_000441.2).
Immunogen Sequence	WSYNTSTEAMTYDEASAYCQ QRYTHLVAIONKEEIEYLNS ILSYSPSYWIGIRKVNNVW VVVGTKPLTEEAKNWAPGE
Sequence	PNNRQKDEDCVEIYIKREKD VGMWNDERCSKKLALCYTA ACTNTSCSGHGECVETINNY TCKCDPGFSGLKCEQIVNCT ALESPEHGSLVCSHPLGNFS YNSSCSISCDRGYLPSSMET MQCMSSGEW