

Anti-LCK antibody (1-210 aa) [SRM] (STJ11107321)

STJ11107321

GENERAL INFORMATION



Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:1000-1:4000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	3932
Gene Symbol	LCK
Uniprot ID	LCK_HUMAN
Immunogen	
Immunogen Region	1-210 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human LCK (NP_005347.3). MGCGCSSHPEDDWMENIDVC ENCHYIVPLDGKGTLLIRN GSEVRDPLVTYEGSNPPASP LQDNLVIALHSYEPSHDGDL
Immunogen Sequence	GFEKGEQLRILEQSGEWWKA QSLTTGQEGFIPNFVAKAN SLEPEPWFFKNLSRKDAERQ LLAPGNTHGSFLIRESESTA GSFSLSVRDFDQNGGEVVKH YKIRNLNNGGFYISPRITFP GLHELVRHYT



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081