

Anti-TAP2 antibody (450-686 aa) [SRM] (STJ11107348)

STJ11107348

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:1000-1:4000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	6891
Gene Symbol	TAP2
Uniprot ID	TAP2_HUMAN
Immunogen	
Immunogen Region	450-686 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 450-686 of human TAP2 (NP_001276972.1). QPNLPSPGTLAPTTLQGVVK FQDVSFAYPNRPDRPVKGL TFTLRPGVETALVGPNGSGK STVAALLONLYQPTGGQVLL DEKIPISQYEHICYLHSQVSV GQEPVLFSGSVRNIAIYGLQ SCEDDKVMAAAQAAHADDFI QEMEHGIYTDVGEKGSQLAA GQKQRLAIARALVRDPRVLI LDEATSALDVQCEQALQDWN SRGDRTVLVIAHRLQTVQRA HQILVLQEGKLQKLAQL
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081