

Anti-Synapsin-1 antibody (401-500 aa) [SRM] (STJ11107349)

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IF/ICC: 1:200-1:800 ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	6853
Gene Symbol	SYN1
Uniprot ID	SYN1_HUMAN
Immunogen	
Immunogen Region	401-500 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 401-500 of human Synapsin-1 (NP_008881.2 ,).
Immunogen Sequence	EDKQLIVLVNKMALPR QRQRDASPRGSHGQTPSPG ALPLGRQTSQQPAGPPAQQR PPPQGGPPQPGPGPQRQGPP LQQRPPPPQGQQLSGLGPPA