

Anti-VIP antibody (31-170 aa) [SRM] (STJ11107351)

STJ11107351

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	IF/ICC:1:200-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	7432
Gene Symbol	VIP
Uniprot ID	VIP_HUMAN
Immunogen	
Immunogen Region	31-170 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 31-170 of human VIP (NP_003372.1).
Immunogen Sequence	APSA LR LGDRIPFEGANEPD QVSLKEDIDMLQNALAENDT PYYDVSRNARHADGVFTSDF SKLLGQLSAKKYLESLMGKR VSSNISED PV PKRHSDAVF TDNYTRLRKQMAVKKYLN SI LNGKRSSEGESPD F PEELEK