

Anti-CD304/Neuropilin-1 antibody (22-856 aa) [SRM] (STJ11107373)
STJ11107373

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:1000-1:2000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	18186
Gene Symbol	Nrp1
Uniprot ID	NRP1_MOUSE
Immunogen	
Immunogen Region	22-856 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 22-856 of mouse CD304/Neuropilin-1 (NP_032763.2).
Immunogen Sequence	FRSDKCGGTIKENPGYLTS PGYPHSYHPSEKCEWLQAP EPYQRIMINFNPHFDLEDRD CKYDYVEIDGENEGGRLWG KFCGKIAPSPVSSGPFLFI KVVSDYETHGAGFSIRYEIF KRGPECSQNYTAPTGVKSP GFPEKYPNSELECTYIIIFAPK MSEIILEFESFDLEQDSNPP GGMFCRYDRLEIWDGFPFVG PHIGRYCGQKTPGRIRSSSG VLSMVFYTDSAIAKEGFSA



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081