

Anti-KCNA5 antibody (1-120 aa) (STJ11107391)

STJ11107391

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	3741
Gene Symbol	KCNA5
Uniprot ID	KCNA5_HUMAN
Immunogen	
Immunogen Region	1-120 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human KCNA5 (NP_002225.2).
Immunogen Sequence	MEIALVPLENGGAMTVRGGD EARAGCGQATGGELQCPPTA GLSDGPKEPAPKGRGAQRDA DSGVRPLPLPDGVRPLPP LPEELPRPRRPPPEDEEEEG DPGLGTVEDQALGTASLHHQ